
edtools

Release 1.0.4

Stef Smeets

Aug 20, 2022

API REFERENCE

1	Installation	3
2	OS Requirement	5
3	Software Requirements	7
4	Package dependencies	9
5	Documentation	11
6	Pipeline tools	13
6.1	autoindex.py	13
6.2	extract_xds_info.py	13
6.3	find_cell.py	14
6.4	make_xscale.py	14
6.5	cluster.py	14
7	Helper tools	15
7.1	make_shelx.py	15
7.2	run_pointless.py	15
7.3	update_xds.py	15
7.4	find_rotation_axis.py	16
8	Demo of using edtools to process batch 3D electron diffraction datasets	17
9	API Reference	19
9.1	edtools	19
10	Examples	41
10.1	edtools Demo	41
10.2	Instruction for using on your own data	59
11	Links	61
12	Indices and tables	63
	Python Module Index	65
	Index	67

Collection of tools for automated processing and clustering of batch 3-dimensional electron diffraction (3D ED) datasets.

The source for this project is available [here](#).

INSTALLATION

Install using `pip install edtools`. Installation should take less than 20 seconds on a normal desktop.

Find the latest [releases](#) for the versions that have been tested on.

OS REQUIREMENT

Windows 10 or newer.

SOFTWARE REQUIREMENTS

- Python 3.6+ including `numpy`, `scipy`, `matplotlib`, and `pandas` libraries
- ``sginfo`` <<https://github.com/rwgk/sginfo>>`__ or ``cctbx.python`` <<https://cctbx.github.io/installation.html#installation>>`__ must be available on the system path for `edtools.make_shelx`
- Access to `WSL`
- XDS package must be installed properly under `WSL`

PACKAGE DEPENDENCIES

Check [pyproject.toml](#) for the full dependency list and versions.

DOCUMENTATION

See the documentation at <https://edtools.readthedocs.io>.

PIPELINE TOOLS

At any step, run `edtools.xxx -h` for help with possible arguments.

6.1 `autoindex.py`

Looks for files matching `XDS.INP` in all subdirectories and runs them using `XDS`.

```
In: XDS.INP
Out: XDS data processing on all files
```

Usage:

```
edtools.autoindex
```

6.2 `extract_xds_info.py`

Looks files matching `CORRECT.LP` in all subdirectories and extracts unit cell/integration info. Summarizes the unit cells in the excel file `cells.xlsx` and `cells.yaml`. `XDS_ASCII.HKL` files matching the completeness / `CC(1/2)` criteria are listed in `filelist.txt`. Optionally, gathers the corresponding `XDS_ASCII.HKL` files in the local directory. The `cells.yaml` file can be used as input for further processing.

```
In: CORRECT.LP
Out: cells.yaml
     cells.xlsx
     filelist.txt
```

Usage:

```
edtools.extract_xds_info
```

6.3 find_cell.py

This program a cells.yaml file and shows histogram plots with the unit cell parameters. This program mimicks `CELLPARM` <http://xds.mpimf-heidelberg.mpg.de/html_doc/cellparm_program.html> and calculates the weighted mean lattice parameters, where the weight is typically the number of observed reflections (defaults to 1.0). For each lattice parameter, the mean is calculated in a given range (default range = median+-2). The range can be changed by dragging the cursor on the histogram plots.

Alternatively, the unit cells can be clustered by giving the `--cluster` command, in which a dendrogram is shown. The cluster cutoff can be selected by clicking in the dendrogram. The clusters will be written to `cells_cluster_#.yaml`.

```
In:  cells.yaml
Out: mean cell parameters
     cells_*.yaml (clustering only)
```

Usage:

```
edtools.find_cell cells.yaml --cluster
```

6.4 make_xscale.py

Prepares an input file `XSCALE.INP` for `XSCALE` and corresponding `XDSCONV.INP` for `XDSCONV`. Takes a `cells.yaml` file or a series of `XDS_ASCII.HKL` files as input, and uses those to generate the `XSCALE.INP` file.

```
In:  cells.yaml / XDS_ASCII.HKL
Out: XSCALE.INP
```

Usage:

```
edtools.make_xscale cells.yaml -c 10.0 20.0 30.0 90.0 90.0 90.0 -s Cmmm
```

6.5 cluster.py

Parses the `XSCALE.LP` file for the correlation coefficients between reflection files to perform hierarchical cluster analysis (Giordano et al., Acta Cryst. (2012). D68, 649–658). The cutoff threshold can be selected by clicking in the dendrogram window. The program will write new `XSCALE.LP` files to subdirectories `cluster_#`, and run `XSCALE` on them, and (if available), pointless.

```
In:  XSCALE.LP
Out: cluster_n/
     filelist.txt
     *_XDS_ASCII.HKL
     XSCALE processing
     Pointless processing
     shelx.hkl
     shelx.ins (optional)
```

Usage:

```
edtools.cluster
```

HELPER TOOLS

7.1 make_shelx.py

Creates a shelx input file. Requires `sginfo` to be available on the system path to generate the SYMM/LATT cards.

In: cell, space group, composition
Out: shelx.ins

Usage:

```
edtools.make_shelx -c 10.0 20.0 30.0 90.0 90.0 90.0 -s Cmmm -m Si180 0360
```

7.2 run_pointless.py

Looks for XDS_ASCII.HKL files specified in the cells.yaml, or on the command line and runs Pointless on them.

In: cells.yaml / XDS_ASCII.HKL
Out: Pointless processing

7.3 update_xds.py

Looks files matching CORRECT.LP in all subdirectories, and updates the cell parameters / space group as specified.

In: XDS.INP
Out: XDS.INP

Usage:

```
edtools.update_xds -c 10.0 20.0 30.0 90.0 90.0 90.0 -s Cmmm
```

7.4 find_rotation_axis.py

Finds the rotation axis and prints out the inputs for several programs (XDS, PETS, DIALS, Instamatic, and RED). Implements the algorithm from Gorelik et al. (Introduction to ADT/ADT3D. In *Uniting Electron Crystallography and Powder Diffraction* (2012), 337-347). The program reads `XDS.INP` to get information about the wavelength, pixelsize, oscillation angle, and beam center, and `SPOT.XDS` (generated by COLSPOT) for the peak positions. If the `XDS.INP` file is not specified, the program will try to look for it in the current directory.

In: <code>XDS.INP</code> , <code>SPOT.XDS</code> Out: Rotation axis
--

Usage:

<code>edtools.find_rotation_axis [XDS.INP]</code>

DEMO OF USING EDTOOLS TO PROCESS BATCH 3D ELECTRON DIFFRACTION DATASETS

See the demo at https://edtools.readthedocs.io/en/latest/examples/edtools_demo.html.

API REFERENCE

9.1 edtools

9.1.1 edtools.autoindex

Functions:

<code>clear_files(path)</code>	Clear LP files
<code>connect(payload)</code>	Try to connect to <i>instamatic</i> indexing server
<code>main()</code>	
<code>parse_xds(path[, sequence])</code>	Parse XDS output (CORRECT.LP) and print summary about indexing progress to the screen.
<code>xds_index(path[, sequence, clear, parallel])</code>	Run XDS at given path.

`edtools.autoindex.clear_files(path: str) → None`
Clear LP files

`edtools.autoindex.connect(payload: str) → None`
Try to connect to *instamatic* indexing server

Parameters

payload (*str*) – Directory where XDS should be run.

`edtools.autoindex.main()`

`edtools.autoindex.parse_xds(path: str, sequence: int = 0) → None`
Parse XDS output (CORRECT.LP) and print summary about indexing progress to the screen.

Parameters

- **path** (*str*) – Path in which XDS has been run
- **sequence** (*int*) – Sequence number, needed for output and house-keeping

`edtools.autoindex.xds_index(path: str, sequence: int = 0, clear: bool = True, parallel: bool = True) → None`
Run XDS at given path.

Parameters

- **path** (*str*) – Run XDS in this directory, expects XDS.INP in this directory
- **sequence** (*int*) – Sequence number, needed for output and house-keeping
- **clear** (*bool*) – Clear some LP files before running XDS

- **parallel** (*bool*) – Call *xds_par* rather than *xds*

9.1.2 edtools.cif_tools

Functions:

add_instrument()

delete_instrument()

list_instrument()

main()

update_cif(*fn*[, *wavelength*, *description*, ...])

`edtools.cif_tools.add_instrument()`

`edtools.cif_tools.delete_instrument()`

`edtools.cif_tools.list_instrument()`

`edtools.cif_tools.main()`

`edtools.cif_tools.update_cif(fn, wavelength=None, description=None, color=None, remove_hkl=None, reply=None, instrument=None)`

9.1.3 edtools.cluster

Functions:

check_for_pointless()

clean_params(*inp*)

distance_from_dendrogram(*z*[, *distance*])

get_clusters(*z*[, *distance*, *fns*, *method*, ...])

get_condensed_distance_matrix(*corrmat*)

main()

parse_xscale_lp(*fn*)

parse_xscale_lp_initial(*[fn]*)

run_pointless(*filepat*[, *verbose*, *i*])

run_xscale(*clusters*, *cell*, *spgr*[, ...])

```

edtools.cluster.check_for_pointless()
edtools.cluster.clean_params(inp)
edtools.cluster.distance_from_dendrogram(z, distance=None)
edtools.cluster.get_clusters(z, distance=0.5, fns=[], method='average', min_size=1)
edtools.cluster.get_condensed_distance_matrix(corrmat)
edtools.cluster.main()
edtools.cluster.parse_xscale_lp(fn)
edtools.cluster.parse_xscale_lp_initial(fn='XSCALE.LP')
edtools.cluster.run_pointless(filepat, verbose=True, i=0)
edtools.cluster.run_xscale(clusters, cell, spgr, resolution=(20.0, 0.8), ioversigma=2)

```

9.1.4 edtools.extract_xds_info

Functions:

<code>cells_to_cellparm(ps)</code>	Takes a list of <code>xds_parser</code> instances and writes the cell parameters to an instruction file <code>CELLPARAM.INP</code> for the program <code>cellparm</code> .
<code>cells_to_excel(ps[, fn])</code>	Takes a list of <code>xds_parser</code> instances and writes the cell parameters to an excel file <code>cells.xlsx</code> .
<code>cells_to_yaml(ps[, fn])</code>	
<code>cells_to_yaml_xparam(uc[, fn])</code>	
<code>evaluate_symmetry(ps)</code>	
<code>gather_xds_ascii(ps[, min_completeness, ...])</code>	Takes a list of <code>xds_parser</code> instances and gathers the corresponding <code>XDS_ASCII.HKL</code> files into the current directory.
<code>lattice_to_space_group(lattice)</code>	
<code>main()</code>	
<code>parse_xparam_for_uc(fn)</code>	

Classes:

<code>xds_parser(filename)</code>	docstring for <code>xds_parser</code>
-----------------------------------	---------------------------------------

```
edtools.extract_xds_info.cells_to_cellparm(ps)
```

Takes a list of `xds_parser` instances and writes the cell parameters to an instruction file `CELLPARAM.INP` for the program `cellparm`.

`edtools.extract_xds_info.cells_to_excel(ps, fn='cells.xlsx')`

Takes a list of *xds_parser* instances and writes the cell parameters to an excel file *cells.xlsx*.

`edtools.extract_xds_info.cells_to_yaml(ps, fn='cells.yaml')`

`edtools.extract_xds_info.cells_to_yaml_xparm(uc, fn='cells_xparm.yaml')`

`edtools.extract_xds_info.evaluate_symmetry(ps)`

`edtools.extract_xds_info.gather_xds_ascii(ps, min_completeness=10.0, min_cchalf=90.0, gather=False)`

Takes a list of *xds_parser* instances and gathers the corresponding *XDS_ASCII.HKL* files into the current directory. The data source and numbering scheme is summarized in the file *filelist.txt*.

`edtools.extract_xds_info.lattice_to_space_group(lattice)`

`edtools.extract_xds_info.main()`

`edtools.extract_xds_info.parse_xparm_for_uc(fn)`

class `edtools.extract_xds_info.xds_parser(filename)`

Bases: object

docstring for *xds_parser*

Methods:

`__init__(filename)`

`cell_as_dict()`

`cell_info([sequence])`

`info_header([hline])`

`integration_info([sequence, outer_shell, ...])`

`parse()`

`print_filename()`

Attributes:

`space_group`

`unit_cell`

`volume`

`__init__(filename)`

`cell_as_dict()`

`cell_info(sequence=0)`

```

static info_header(hline=True)

integration_info(sequence=0, outer_shell=True, filename=False)

parse()

print_filename()

property space_group

property unit_cell

property volume

```

9.1.5 edtools.find_beam_center

Functions:

<code>find_beam_center(img[, sigma, m, kind])</code>	Find the center of the primary beam in the image <i>img</i> . The position is determined by summing along X/Y directions and finding the position along the two directions independently.
<code>find_peak_max(arr, sigma[, m, w, kind])</code>	Find the index of the pixel corresponding to peak maximum in 1D pattern <i>arr</i> .
<code>main()</code>	
<code>read_adsc(fname)</code>	read in the file.
<code>readheader(infile)</code>	read an adsc header.
<code>swap_needed(header)</code>	
<code>translate_image(arr, shift)</code>	Translate an image according to shift.
<code>write_adsc(fname, data[, header])</code>	Write adsc format.

`edtools.find_beam_center.find_beam_center`(*img*: ~numpy.ndarray, *sigma*: int = 30, *m*: int = 100, *kind*: int = 3) -> (<class 'float'>, <class 'float'>)

Find the center of the primary beam in the image *img*. The position is determined by summing along X/Y directions and finding the position along the two directions independently.

Uses interpolation by factor *m* to find the coordinates of the primary beam with subpixel accuracy.

`edtools.find_beam_center.find_peak_max`(*arr*: ~numpy.ndarray, *sigma*: int, *m*: int = 50, *w*: int = 10, *kind*: int = 3) -> (<class 'float'>, <class 'float'>)

Find the index of the pixel corresponding to peak maximum in 1D pattern *arr*.

First, the pattern is smoothed using a gaussian filter with standard deviation *sigma*. The initial guess takes the position corresponding to the largest value in the resulting pattern. A window of size $2*w+1$ around this guess is taken and expanded by factor *m* to interpolate the pattern to get the peak maximum position with subpixel precision.

`edtools.find_beam_center.main`()

`edtools.find_beam_center.read_adsc`(*fname*: str) -> (<built-in function array>, <class 'dict'>)
read in the file.

`edtools.find_beam_center.readheader(infile)`

read an adsc header.

`edtools.find_beam_center.swap_needed(header: dict) → bool`

`edtools.find_beam_center.translate_image(arr, shift: array) → array`

Translate an image according to shift. Shift should be a 2D numpy array

`edtools.find_beam_center.write_adsc(fname: str, data: array, header: dict = {})`

Write adsc format.

9.1.6 edtools.find_cell

Functions:

<code>cluster_cell(cells[, distance, method, ...])</code>	Perform hierarchical cluster analysis on a list of cells.
<code>d_calculator(cell)</code>	Helper function for <code>unit_cell_lcv_distance</code>
<code>distance_from_dendrogram(z[, ylabel, ...])</code>	Takes a linkage object <code>z</code> from <code>scipy.cluster.hierarchy.linkage</code> and displays a dendrogram.
<code>find_cell(cells, weights[, binsize])</code>	Opens a plot with 6 subplots in which the cell parameter histogram is displayed.
<code>get_clusters(z, cells[, distance])</code>	
<code>main()</code>	
<code>parse_cellparm(fn)</code>	
<code>put_in_order(cells)</code>	order cell parameters in order to eliminate difference in cell distance because of parameter order
<code>to_radian(cells)</code>	convert all angles in unit cell parameter list to radians cells: the cell parameters that are parsed from <code>cells.yaml</code> as np array
<code>to_sin(cells)</code>	convert all angles in unit cell parameter list to sine cells: the cell parameters that are parsed from <code>cells.yaml</code> as np array
<code>unit_cell_lcv_distance(cell1, cell2)</code>	Implements Linear Cell Volume from Acta Cryst.
<code>volume_difference(cell1, cell2)</code>	Return the absolute difference in volumes between two unit cells
<code>weighted_average(values[, weights])</code>	Returns weighted mean and standard deviation

`edtools.find_cell.cluster_cell(cells: list, distance: Optional[float] = None, method: str = 'average', metric: str = 'euclidean', use_radian: bool = False, use_sine: bool = False)`

Perform hierarchical cluster analysis on a list of cells.

method: lcv, volume, euclidean distance: cutoff distance, if it is not given, pop up a dendrogram to

interactively choose a cutoff distance

use_radian: Use radian instead of degrees to downweight difference use_sine: Use sine for unit cell clustering (to disambiguate the difference in angles)

`edtools.find_cell.d_calculator(cell: list) → tuple`

Helper function for `unit_cell_lcv_distance`

`edtools.find_cell.distance_from_dendrogram`(*z*, *ylabel*: *str* = "", *initial_distance*: *Optional[float]* = *None*)
→ float

Takes a linkage object *z* from `scipy.cluster.hierarchy.linkage` and displays a dendrogram. The cutoff distance can be picked interactively, and is returned *ylabel*: sets the label for the y-axis *initial_distance*: initial cutoff distance to display

`edtools.find_cell.find_cell`(*cells*, *weights*, *binsize*=0.5)

Opens a plot with 6 subplots in which the cell parameter histogram is displayed. It will calculate the weighted mean of the unit cell parameters. The ranges can be adjusted by dragging on the plots.

`edtools.find_cell.get_clusters`(*z*, *cells*, *distance*=0.5)

`edtools.find_cell.main`()

`edtools.find_cell.parse_cellparm`(*fn*)

`edtools.find_cell.put_in_order`(*cells*)

order cell parameters in order to eliminate difference in cell distance because of parameter order

`edtools.find_cell.to_radian`(*cells*)

convert all angles in unit cell parameter list to radians *cells*: the cell parameters that are parsed from `cells.yaml` as np array

`edtools.find_cell.to_sin`(*cells*)

convert all angles in unit cell parameter list to sine cells: the cell parameters that are parsed from `cells.yaml` as np array

`edtools.find_cell.unit_cell_lcv_distance`(*cell1*: *list*, *cell2*: *list*) → float

Implements Linear Cell Volume from Acta Cryst. (2013). D69, 1617-1632

`edtools.find_cell.volume_difference`(*cell1*: *list*, *cell2*: *list*)

Return the absolute difference in volumes between two unit cells

`edtools.find_cell.weighted_average`(*values*, *weights*=*None*)

Returns weighted mean and standard deviation

9.1.7 edtools.find_rotation_axis

Functions:

<code>cylinder_histo(xyz[, bins])</code>	Take reciprocal lattice vectors in XYZ format and output cylindrical projection.
<code>load_spot_xds(fn, beam_center, osc_angle, ...)</code>	Load the given SPOT.XDS file (<i>fn</i>) and return an array with the reciprocal
<code>main()</code>	
<code>make(arr, omega, wavelength)</code>	Prepare xyz (reciprocal space coordinates) from reflection positions/angle (<i>arr</i>), which is the list of reflections read from XDS (SPOT.XDS)
<code>make_2d_rotmat(theta)</code>	Take angle in radians, and return 2D rotation matrix
<code>optimize(arr, omega_start[, wavelength, ...])</code>	Optimize the value of omega around the given point.
<code>parse_xds_inp(fn)</code>	Parse the XDS.INP file to find the required numbers for the optimization Looks for wavelength, pixelsize, beam_center, oscillation range
<code>plot_histo(H, xedges, yedges[, title])</code>	Plot the histogram of the cylindrical projection.
<code>random_sample(arr, n)</code>	Select random sample of <i>n</i> rows from array
<code>rotation_axis_to_xyz(rotation_axis[, ...])</code>	Convert rotation axis angle to XYZ vector compatible with 'xds', or 'dials' Set invert to 'True' for anti-clockwise rotation
<code>rotation_matrix(axis, theta)</code>	Calculates the rotation matrix around axis of angle theta (radians)
<code>xyz2cyl(arr)</code>	Take a set of reflections in XYZ and convert to polar (cylindrical) coordinates

`edtools.find_rotation_axis.cylinder_histo(xyz, bins=(1000, 500))`

Take reciprocal lattice vectors in XYZ format and output cylindrical projection. *Bins* gives the resolution of the 2D histogram.

`edtools.find_rotation_axis.load_spot_xds(fn, beam_center: [<class 'float'>, <class 'float'>], osc_angle: float, pixelsize: float)`

Load the given SPOT.XDS file (*fn*) and return an array with the reciprocal

x, *y*, and angle for the centroid of each reflection

beam_center: coordinates of the primary beam, read from XDS.INP *osc_angle*: oscillation_angle (degrees) per frame, will be multiplied by the average frame number

that a reflection appears on (column 3 in *arr*)

pixelsize: defined in px/Ångström

http://xds.mpimf-heidelberg.mpg.de/html_doc/xds_files.html#SPOT.XDS

`edtools.find_rotation_axis.main()`

`edtools.find_rotation_axis.make(arr, omega: float, wavelength: float)`

Prepare xyz (reciprocal space coordinates) from reflection positions/angle (*arr*), which is the list of reflections read from XDS (SPOT.XDS)

omega: rotation axis (degrees), which is defined by the angle between x

(horizontal axis pointing right) and the rotation axis going in clockwise direction

Note that:

1. $x \leftrightarrow y$ are flipped

This is to ensure to match the XDS convention with the one I'm used to

`edtools.find_rotation_axis.make_2d_rotmat(theta)`

Take angle in radians, and return 2D rotation matrix

`edtools.find_rotation_axis.optimize(arr, omega_start: float, wavelength=<class 'float'>, plusminus: int = 180, step: int = 10, hist_bins: (<class 'int'>, <class 'int'>) = (1000, 500), plot: bool = False) → float`

Optimize the value of omega around the given point.

`omega_start`: defines the starting angle step, `plusminus`: together with `omega_start` define the range of values to loop over `hist_bins`: size of the 2d histogram to produce the final phi/theta plot `plot`: toggle to plot the histogram after each step

`edtools.find_rotation_axis.parse_xds_inp(fn)`

Parse the XDS.INP file to find the required numbers for the optimization Looks for wavelength, pixelsize, beam_center, oscillation range

`edtools.find_rotation_axis.plot_histo(H, xedges, yedges, title='Histogram')`

Plot the histogram of the cylindrical projection.

`edtools.find_rotation_axis.random_sample(arr, n)`

Select random sample of n rows from array

`edtools.find_rotation_axis.rotation_axis_to_xyz(rotation_axis, invert=False, setting='xds')`

Convert rotation axis angle to XYZ vector compatible with 'xds', or 'dials' Set `invert` to 'True' for anti-clockwise rotation

`edtools.find_rotation_axis.rotation_matrix(axis, theta)`

Calculates the rotation matrix around axis of angle theta (radians)

`edtools.find_rotation_axis.xyz2cyl(arr)`

Take a set of reflections in XYZ and convert to polar (cylindrical) coordinates

9.1.8 edtools.make_shelx

9.1.9 edtools.make_xscale

Functions:

`get_xds_ascii_names(lst)`

`main()`

`parse_xds_ascii(fn)`

`write_xdsconv_inp(resolution)`

`write_xscale_inp(fns, unit_cell, ...)`

`edtools.make_xscale.get_xds_ascii_names(lst)`

`edtools.make_xscale.main()`

`edtools.make_xscale.parse_xds_ascii(fn)`

`edtools.make_xscale.write_xdsconv_inp(resolution)`

`edtools.make_xscale.write_xscale_inp(fns, unit_cell, space_group, resolution)`

9.1.10 edtools.period_table

Classes:

`App(parent, *args, **kwargs)`

Functions:

`main()`

class `edtools.period_table.App(parent, *args, **kwargs)`

Bases: `Frame`

Methods:

<code>__init__(parent, *args, **kwargs)</code>	Construct a frame widget with the parent MASTER.
<code>info(text)</code>	Return information about the packing options for this widget.
<code>name(text)</code>	

`__init__(parent, *args, **kwargs)`

Construct a frame widget with the parent MASTER.

Valid resource names: background, bd, bg, borderwidth, class, colormap, container, cursor, height, highlightbackground, highlightcolor, highlightthickness, relief, takefocus, visual, width.

info(*text*)

Return information about the packing options for this widget.

name(*text*)

`edtools.period_table.main()`

9.1.11 edtools.reflection_tool

Classes:

<i>GroupReflectionsGUI</i> (parent)	A GUI frame for reflections grouping
-------------------------------------	--------------------------------------

Functions:

main()

class edtools.reflection_tool.**GroupReflectionsGUI**(*parent*)

Bases: Labelframe

A GUI frame for reflections grouping

Methods:

<code>__init__(parent)</code>	Construct a Tk Labelframe with parent master.
<code>calc_structure_factors(structures[, dmin, ...])</code>	Wrapper around <code>f_calc_structure_factors()</code> Takes a structure object in which there is only one structure
<code>check_I_frame_seq()</code>	
<code>corr_prec()</code>	
<code>exti_corr(value, power, param)</code>	
<code>exti_corr_der(value, power, param)</code>	
<code>f_calc_structure_factors(structure, **kwargs)</code>	Takes cctbx structure and returns f_calc miller array Takes an optional options dictionary with keys: input: **kwargs : 'd_min': minimum d-spacing for structure factor calculation 'algorithm': which algorithm to use ('direct', 'fft', 'automatic') structure: <cctbx.xray.structure.structure object> output: f_calc: <cctbx.miller.array object> with calculated structure factors in the <code>f_calc.data()</code> function
<code>gen_fcalc()</code>	
<code>group_df()</code>	
<code>init_vars()</code>	
<code>lorentz_corr(df)</code>	
<code>open_file()</code>	
<code>point_to_rotation_axis(x, y, center, slope)</code>	
<code>read_cif(f)</code>	
<code>remove_reflection()</code>	
<code>save_file(df)</code>	
<code>save_grouped()</code>	
<code>scaling_factor(value, power)</code>	
<code>scaling_func(value, power)</code>	
<code>split_grouped()</code>	
<code>transform_integrated()</code>	
<code>validate(action, index, value_if_allowed, ...)</code>	
<code>validate_range(action, index, ...)</code>	
<code>__init__(parent)</code>	

Construct a Ttk Labelframe with parent master.

STANDARD OPTIONS

class, cursor, style, takefocus

WIDGET-SPECIFIC OPTIONS

labelanchor, text, underline, padding, labelwidget, width, height

calc_structure_factors(*structures*, *dmin=1.0*, *table='electron'*, *prefix=""*, *verbose=True*, ***kwargs*)

Wrapper around f_calc_structure_factors() Takes a structure object in which there is only one structure

dmin can be a dataframe and it will take the minimum dspacing (as specified by col 'd') or a float if combine is specified, function will return a dataframe combined with the given one, otherwise a dictionary of dataframes

prefix is a prefix for the default names fcalc/phases to identify different structures

check_I_frame_seq()

corr_prec()

exti_corr(*value*, *power*, *param*)

exti_corr_der(*value*, *power*, *param*)

f_calc_structure_factors(*structure*, ***kwargs*)

Takes cctbx structure and returns f_calc miller array Takes an optional options dictionary with keys: input:

****kwargs:**

'd_min': minimum d-spacing for structure factor calculation 'algorithm': which algorithm to use ('direct', 'fft', 'automatic')

structure: <cctbx.xray.structure.structure object>

output:

f_calc: <cctbx.miller.array object> with calculated structure factors
in the f_calc.data() function

gen_fcalc()

group_df()

init_vars()

lorentz_corr(*df*)

open_file()

point_to_rotation_axis(*x*, *y*, *center*, *slope*)

read_cif(*f*)

remove_reflection()

save_file(*df*)

save_grouped()

`scaling_factor(value, power)`

`scaling_func(value, power)`

`split_grouped()`

`transform_integrated()`

`validate(action, index, value_if_allowed, prior_value, text, validation_type, trigger_type, widget_name)`

`validate_range(action, index, value_if_allowed, prior_value, text, validation_type, trigger_type, widget_name)`

`edtools.reflection_tool.main()`

9.1.12 edtools.run_pointless

Functions:

`main()`

`edtools.run_pointless.main()`

9.1.13 edtools.scattering_factor

Classes:

`FitResult_4p(a0, b0, a1, b1, a2, b2, a3, b3, c)`

`FitResult_5p(a0, b0, a1, b1, a2, b2, a3, b3, ...)`

`ScatteringFactorGUI(parent)` A GUI frame for scattering factor

Functions:

`fit_4_param(s, target[, is_xray, method, ...])`

`fit_5_param(s, target[, is_xray, method, ...])`

`main()`

`class edtools.scattering_factor.FitResult_4p(a0, b0, a1, b1, a2, b2, a3, b3, c)`

Bases: `tuple`

Attributes:

<i>a0</i>	Alias for field number 0
<i>a1</i>	Alias for field number 2
<i>a2</i>	Alias for field number 4
<i>a3</i>	Alias for field number 6
<i>b0</i>	Alias for field number 1
<i>b1</i>	Alias for field number 3
<i>b2</i>	Alias for field number 5
<i>b3</i>	Alias for field number 7
<i>c</i>	Alias for field number 8

property a0

Alias for field number 0

property a1

Alias for field number 2

property a2

Alias for field number 4

property a3

Alias for field number 6

property b0

Alias for field number 1

property b1

Alias for field number 3

property b2

Alias for field number 5

property b3

Alias for field number 7

property c

Alias for field number 8

```
class edtools.scattering_factor.FitResult_5p(a0, b0, a1, b1, a2, b2, a3, b3, a4, b4, c)
```

Bases: tuple

Attributes:

<i>a0</i>	Alias for field number 0
<i>a1</i>	Alias for field number 2
<i>a2</i>	Alias for field number 4
<i>a3</i>	Alias for field number 6
<i>a4</i>	Alias for field number 8
<i>b0</i>	Alias for field number 1
<i>b1</i>	Alias for field number 3
<i>b2</i>	Alias for field number 5
<i>b3</i>	Alias for field number 7
<i>b4</i>	Alias for field number 9
<i>c</i>	Alias for field number 10

property a0

Alias for field number 0

property a1

Alias for field number 2

property a2

Alias for field number 4

property a3

Alias for field number 6

property a4

Alias for field number 8

property b0

Alias for field number 1

property b1

Alias for field number 3

property b2

Alias for field number 5

property b3

Alias for field number 7

property b4

Alias for field number 9

property c

Alias for field number 10

class edtools.scattering_factor.ScatteringFactorGUI(*parent*)

Bases: Labelframe

A GUI frame for scattering factor

Methods:

<code>__init__(parent)</code>	Construct a Tk Labelframe with parent master.
<code>clear_plot()</code>	
<code>draw()</code>	
<code>fit_4_param()</code>	
<code>fit_5_param()</code>	
<code>func_4p(param, s)</code>	
<code>func_4p_charge(param, s[, charge])</code>	
<code>func_5p(param, s)</code>	
<code>func_5p_charge(param, s[, charge])</code>	
<code>get_parameters(scat_lib, source, atom, charge)</code>	
<code>init_vars()</code>	
<code>load_lib()</code>	
<code>output()</code>	
<code>plot()</code>	
<code>read()</code>	
<code>select_element(element)</code>	

__init__(parent)

Construct a Tk Labelframe with parent master.

STANDARD OPTIONS

class, cursor, style, takefocus

WIDGET-SPECIFIC OPTIONS

labelanchor, text, underline, padding, labelwidget, width, height

clear_plot()**draw()****fit_4_param()****fit_5_param()****func_4p(param, s)****func_4p_charge(param, s, charge=0)**

`func_5p(param, s)`
`func_5p_charge(param, s, charge=0)`
`get_parameters(scats_lib, source, atom, charge)`
`init_vars()`
`load_lib()`
`output()`
`plot()`
`read()`
`select_element(element)`

`edtools.scattering_factor.fit_4_param(s, target, is_xray=False, method: str = 'leastsq', verbose: bool = False, **param)`

`edtools.scattering_factor.fit_5_param(s, target, is_xray=False, method: str = 'leastsq', verbose: bool = False, **param)`

`edtools.scattering_factor.main()`

9.1.14 edtools.update_xds

Functions:

`main()`

`update_xds(fn[, cell, spgr, comment, ...])`

`edtools.update_xds.main()`

`edtools.update_xds.update_xds(fn, cell=None, spgr=None, comment=False, axis_error=None, angle_error=None, overload=None, lo_res=None, hi_res=None, cut_frames=None, wfacI=None, apd=None, jobs=None, sp=None, indnumthre=None, d=False, dl=None, processors=None, center=None, axis=None, cam_len=None, mosaicity=None, pixel_size=None, untrusted=None, corr=None, refine_idx=None, refine_integrate=None, refine_corr=None, trusted_region=None, trusted_pixels=None)`

9.1.15 edtools.utils

Functions:

<code>parse_args_for_fns(args[, name, match])</code>	Parse list of filenames and resolve wildcards name: Name of the file to locate match: Match the file list against the provided glob-style pattern. If the match is False, the path is removed from the list. example: match="SMV_reprocessed".
<code>space_group_lib()</code>	Initialize simple space group library mapping the space group number to a dict with information on the <i>class</i> (crystal class), <i>lattice</i> (lattice symbol), <i>laue_symmetry</i> (number of the lowest symmetry space group for this lattice), <i>name</i> (space group name), and <i>number</i> (space group number).
<code>volume(cell)</code>	Returns volume for the general case from cell parameters

`edtools.utils.parse_args_for_fns(args, name='XDS.INP', match=None)`

Parse list of filenames and resolve wildcards name:

Name of the file to locate

match:

Match the file list against the provided glob-style pattern. If the match is False, the path is removed from the list. example:

match="SMV_reprocessed"

`edtools.utils.space_group_lib()`

Initialize simple space group library mapping the space group number to a dict with information on the *class* (crystal class), *lattice* (lattice symbol), *laue_symmetry* (number of the lowest symmetry space group for this lattice), *name* (space group name), and *number* (space group number).

`edtools.utils.volume(cell)`

Returns volume for the general case from cell parameters

9.1.16 edtools.widgets

Classes:

<code>Hoverbox(widget[, text])</code>	create a tooltip for a given widget
<code>Spinbox([master])</code>	Ttk Spinbox is an Entry with increment and decrement arrows It is commonly used for number entry or to select from a list of string values.

class `edtools.widgets.Hoverbox(widget, text='widget info')`

Bases: object

create a tooltip for a given widget

Methods:

`__init__(widget[, text])`

`enter([event])`

`hidetip()`

`leave([event])`

`schedule()`

`showtip([event])`

`unschedule()`

`__init__(widget, text='widget info')``enter(event=None)``hidetip()``leave(event=None)``schedule()``showtip(event=None)``unschedule()`**class** edtools.widgets.**Spinbox**(master=None, **kw)

Bases: Entry

Ttk Spinbox is an Entry with increment and decrement arrows. It is commonly used for number entry or to select from a list of string values.

Methods:

<code>__init__([master])</code>	Construct a Ttk Spinbox widget with the parent master.
---------------------------------	--

<code>set(value)</code>	Sets the value of the Spinbox to value.
-------------------------	---

`__init__(master=None, **kw)`

Construct a Ttk Spinbox widget with the parent master.

STANDARD OPTIONS: class, cursor, style, takefocus, validate, validatecommand, xscrollcommand, invalidcommand

WIDGET-SPECIFIC OPTIONS: to, **from_**, increment, values, wrap, format, command`set(value)`

Sets the value of the Spinbox to value.

9.1.17 edtools.wsl

EXAMPLES

10.1 *edtools* Demo

edtools is a python package for automated processing of a large number of 3D electron diffraction (3D ED) datasets. It can be downloaded from <https://doi.org/10.5281/zenodo.6952810>.

For running *edtools*, *XDS* package for reduction of 3D ED datasets is required. *XDS* package is available at https://xds.mr.mpg.de/html_doc/downloading.html.

A typical cycle of using *edtools* for processing batch 3D ED datasets goes through the following steps:

- `edtools.autoindex`
- `edtools.extract_xds_info`
- `edtools.find_cell`
- `edtools.update_xds`
- `edtools.make_xscale`
- `edtools.cluster`

Here we demonstrate the processing of batch 3D ED datasets for phase analysis and structure determination using *edtools*. The datasets for the demo can be downloaded from https://zenodo.org/record/6533426#.YnoQ7_hBxaQ.

The datasets were collected on a zeolite mixture sample using serial rotation electron diffraction (SerialRED) data collection technique implemented in the program **Instamatic** (available at <https://doi.org/10.5281/zenodo.5175957>), which runs on a JEOL JEM-2100-LaB6 at 200 kV equipped with a 512 x 512 Timepix hybrid pixel detector (55 x 55 μm pixel size, QTPX-262k, Amsterdam Scientific Instruments).

The zeolite mixture sample contains phases **IWVRTH**, and ***CTH**. The information of these three phases can be found from the structure database of zeolites (https://europe.iza-structure.org/IZA-SC/ftc_table.php).

This demo takes around 5-10 min to run on a normal desktop computer with all the required packages installed properly beforehand.

10.1.1 Indexing

Automatically index the 3D ED datasets by running *XDS* in all subfolders (SMV) that contains file *XDS.INP*, which is automatically generated during data collection using *Instamatic*.

[1]: !edtools.autoindex

```
!!! ERROR !!! WRONG TYPE OF INPUT FILE SPECIFIED
!!! ERROR !!! WRONG TYPE OF INPUT FILE SPECIFIED
16 files named XDS.INP (subdir: None) found.
```

```
0: C:\demo\edtools_demo_data\stagepos_0067\crystal_0001\SMV # Mon Aug 1 21:00:56
↪2022
```

```
Spgr 5 - Cell 26.93 14.05 5.36 90.00 90.89 90.00 - Vol ↪
↪2027.80
```

#	dmax	dmin	ntot	nuniq	compl	i/sig	rmeas	CC(1/2)	ISa	B(ov)
0	4.35	0.80	583	324	15.0	4.59	13.7	99.0	7.47	6.72
-	0.85	0.80	54	42	12.5	1.96	26.8	91.9		

```
1: C:\demo\edtools_demo_data\stagepos_0164\crystal_0000\SMV # Mon Aug 1 21:00:58
↪2022
```

```
Spgr 1 - Cell 9.49 9.90 12.47 66.56 89.45 86.35 - Vol ↪
↪1072.59
```

#	dmax	dmin	ntot	nuniq	compl	i/sig	rmeas	CC(1/2)	ISa	B(ov)
1	6.39	0.80	229	209	4.8	24.01	22.3	96.3	50.00	4.74
-	0.91	0.85	31	29	4.5	12.16	21.3	0.0		

```
3: C:\demo\edtools_demo_data\stagepos_0299\crystal_0001\SMV # Mon Aug 1 21:01:00
↪2022
```

```
Spgr 1 - Cell 4.83 14.83 16.03 115.66 89.61 94.16 - Vol ↪
↪1031.87
```

#	dmax	dmin	ntot	nuniq	compl	i/sig	rmeas	CC(1/2)	ISa	B(ov)
3	2.05	0.80	400	312	7.5	2.44	20.6	95.2	4.24	6.11
-	0.84	0.80	27	26	4.7	1.79	11.2	0.0		

```
4: C:\demo\edtools_demo_data\stagepos_0325\crystal_0000\SMV # Mon Aug 1 21:01:02
↪2022
```

```
Spgr 5 - Cell 13.69 25.42 14.90 90.00 115.84 90.00 - Vol ↪
↪4666.87
```

#	dmax	dmin	ntot	nuniq	compl	i/sig	rmeas	CC(1/2)	ISa	B(ov)
4	11.09	0.79	3744	2147	42.8	3.44	13.1	99.6	13.90	8.10
-	0.97	0.90	623	336	47.9	1.32	68.9	84.7		

```
5: C:\demo\edtools_demo_data\stagepos_0341\crystal_0000\SMV # Mon Aug 1 21:01:03
↪2022
```

(continues on next page)

(continued from previous page)

↪2022
 Spgr 5 - Cell 25.67 13.50 17.73 90.00 132.44 90.00 - Vol ↵
 ↪4534.43

#	dmax	dmin	ntot	nuniq	compl	i/sig	rmeas	CC(1/2)	ISa	B(ov)
5	6.86	0.80	2161	1081	21.8	4.23	10.5	99.9	33.30	8.74
-	0.97	0.90	342	159	23.3	0.83	130.3	69.6		

6: C:\demo\edtools_demo_data\stagepos_0368\crystal_0001\SMV # Mon Aug 1 21:01:05↵

↪2022
 Spgr 1 - Cell 10.17 10.36 12.16 93.71 113.40 98.01 - Vol ↵
 ↪1154.16

#	dmax	dmin	ntot	nuniq	compl	i/sig	rmeas	CC(1/2)	ISa	B(ov)
6	10.17	0.80	611	443	9.4	3.17	14.9	97.5	5.07	4.64
-	0.85	0.80	56	53	7.0	1.96	73.1	6.4		

7: C:\demo\edtools_demo_data\stagepos_0538\crystal_0000\SMV # Mon Aug 1 21:01:06↵

↪2022
 Spgr 1 - Cell 10.55 10.52 11.81 80.39 66.60 75.74 - Vol ↵
 ↪1162.33

#	dmax	dmin	ntot	nuniq	compl	i/sig	rmeas	CC(1/2)	ISa	B(ov)
7	5.10	0.80	443	330	7.0	3.80	10.7	99.4	8.61	5.62
-	0.85	0.80	38	36	4.8	1.80	76.5	0.0		

8: C:\demo\edtools_demo_data\stagepos_0648\crystal_0001\SMV # Mon Aug 1 21:01:07↵

↪2022
 Spgr 1 - Cell 13.82 14.32 16.18 86.20 111.75 116.39 - Vol ↵
 ↪2645.41

#	dmax	dmin	ntot	nuniq	compl	i/sig	rmeas	CC(1/2)	ISa	B(ov)
8	6.37	0.80	1460	989	9.1	2.88	16.3	97.8	5.24	7.62
-	0.85	0.80	166	125	7.3	1.36	62.8	56.0		

9: C:\demo\edtools_demo_data\stagepos_0849\crystal_0000\SMV # Mon Aug 1 21:01:08↵

↪2022
 Spgr 5 - Cell 15.06 26.22 15.41 90.00 118.30 90.00 - Vol ↵
 ↪5357.50

#	dmax	dmin	ntot	nuniq	compl	i/sig	rmeas	CC(1/2)	ISa	B(ov)
9	13.11	0.79	2063	1319	22.1	3.46	10.5	99.6	12.09	7.58
-	0.89	0.83	326	223	24.5	1.01	53.6	83.4		

10: C:\demo\edtools_demo_data\stagepos_0905\crystal_0000\SMV # Mon Aug 1 21:01:10↵

↪2022
 Spgr 3 - Cell 13.91 5.07 14.97 90.00 117.96 90.00 - Vol ↵

(continues on next page)

(continued from previous page)

↪932.53

#	dmax	dmin	ntot	nuniq	compl	i/sig	rmeas	CC(1/2)	ISa	B(ov)
10	12.33	0.80	479	300	13.8	3.68	13.4	99.6	16.07	9.46
-	1.20	1.07	58	35	14.6	4.72	24.9	88.5		

11: C:\demo\edtools_demo_data\stagepos_0905\crystal_0001\SMV # Mon Aug 1 21:01:11

↪2022

Spgr 1 - Cell 13.71 14.57 15.77 83.07 68.29 62.34 - Vol ↪
 ↪2587.36

#	dmax	dmin	ntot	nuniq	compl	i/sig	rmeas	CC(1/2)	ISa	B(ov)
11	11.49	0.80	1596	1144	10.7	3.30	12.4	98.5	7.24	7.18
-	0.85	0.80	124	121	7.0	0.94	22.6	83.4		

12: C:\demo\edtools_demo_data\stagepos_0980\crystal_0000\SMV # Mon Aug 1 21:01:13

↪2022

Spgr 1 - Cell 14.56 15.00 15.27 97.22 105.97 120.36 - Vol ↪
 ↪2621.77

#	dmax	dmin	ntot	nuniq	compl	i/sig	rmeas	CC(1/2)	ISa	B(ov)
12	7.54	0.80	1746	1222	11.3	4.00	13.3	98.6	8.77	5.85
-	0.85	0.80	164	146	8.4	1.48	36.5	86.5		

13: C:\demo\edtools_demo_data\stagepos_1014\crystal_0000\SMV # Mon Aug 1 21:01:14

↪2022

Spgr 1 - Cell 5.30 14.56 15.04 112.06 93.44 86.65 - Vol ↪
 ↪1072.87

#	dmax	dmin	ntot	nuniq	compl	i/sig	rmeas	CC(1/2)	ISa	B(ov)
13	5.01	0.81	447	328	7.5	4.11	10.9	98.5	6.65	6.67
-	0.85	0.80	51	44	6.3	2.11	18.7	92.7		

15: C:\demo\edtools_demo_data\stagepos_1283\crystal_0001\SMV # Mon Aug 1 21:01:17

↪2022

Spgr 1 - Cell 13.64 15.02 25.09 93.07 91.13 114.33 - Vol ↪
 ↪4672.25

#	dmax	dmin	ntot	nuniq	compl	i/sig	rmeas	CC(1/2)	ISa	B(ov)
15	6.60	0.80	3124	2149	11.3	3.54	8.4	99.5	12.64	6.94
-	0.85	0.80	346	280	9.2	1.24	56.2	84.9		

10.1.2 Extract cell

Extract the determined unit cell parameters from the output files (CORRECT.LP) of XDS

```
[2]: !edtools.extract_xds_info
14 files named CORRECT.LP (subdir: None) found.
  1: C:\demo\edtools_demo_data\stagepos_0067\crystal_0001\SMV # Mon Aug  1 21:00:56
↪2022
Spgr  5 - Cell      26.93   14.05    5.36   90.00   90.89   90.00 - Vol  ↪
↪2027.80

  2: C:\demo\edtools_demo_data\stagepos_0164\crystal_0000\SMV # Mon Aug  1 21:00:58
↪2022
Spgr  1 - Cell      9.49    9.90   12.47   66.56   89.45   86.35 - Vol  ↪
↪1072.59

  3: C:\demo\edtools_demo_data\stagepos_0299\crystal_0001\SMV # Mon Aug  1 21:01:00
↪2022
Spgr  1 - Cell      4.83   14.83   16.03  115.66   89.61   94.16 - Vol  ↪
↪1031.87

  4: C:\demo\edtools_demo_data\stagepos_0325\crystal_0000\SMV # Mon Aug  1 21:01:02
↪2022
Spgr  5 - Cell     13.69   25.42   14.90   90.00  115.84   90.00 - Vol  ↪
↪4666.87

  5: C:\demo\edtools_demo_data\stagepos_0341\crystal_0000\SMV # Mon Aug  1 21:01:03
↪2022
Spgr  5 - Cell     25.67   13.50   17.73   90.00  132.44   90.00 - Vol  ↪
↪4534.43

  6: C:\demo\edtools_demo_data\stagepos_0368\crystal_0001\SMV # Mon Aug  1 21:01:05
↪2022
Spgr  1 - Cell     10.17   10.36   12.16   93.71  113.40   98.01 - Vol  ↪
↪1154.16

  7: C:\demo\edtools_demo_data\stagepos_0538\crystal_0000\SMV # Mon Aug  1 21:01:06
↪2022
Spgr  1 - Cell     10.55   10.52   11.81   80.39   66.60   75.74 - Vol  ↪
↪1162.33

  8: C:\demo\edtools_demo_data\stagepos_0648\crystal_0001\SMV # Mon Aug  1 21:01:07
↪2022
Spgr  1 - Cell     13.82   14.32   16.18   86.20  111.75  116.39 - Vol  ↪
↪2645.41

  9: C:\demo\edtools_demo_data\stagepos_0849\crystal_0000\SMV # Mon Aug  1 21:01:08
↪2022
Spgr  5 - Cell     15.06   26.22   15.41   90.00  118.30   90.00 - Vol  ↪
↪5357.50

 10: C:\demo\edtools_demo_data\stagepos_0905\crystal_0000\SMV # Mon Aug  1 21:01:10
↪2022
```

(continues on next page)

(continued from previous page)

Spgr 3 - Cell 13.91 5.07 14.97 90.00 117.96 90.00 - Vol 1
 ↪932.53

11: C:\demo\edtools_demo_data\stagepos_0905\crystal_0001\SMV # Mon Aug 1 21:01:11
 ↪2022

Spgr 1 - Cell 13.71 14.57 15.77 83.07 68.29 62.34 - Vol 1
 ↪2587.36

12: C:\demo\edtools_demo_data\stagepos_0980\crystal_0000\SMV # Mon Aug 1 21:01:13
 ↪2022

Spgr 1 - Cell 14.56 15.00 15.27 97.22 105.97 120.36 - Vol 1
 ↪2621.77

13: C:\demo\edtools_demo_data\stagepos_1014\crystal_0000\SMV # Mon Aug 1 21:01:14
 ↪2022

Spgr 1 - Cell 5.30 14.56 15.04 112.06 93.44 86.65 - Vol 1
 ↪1072.87

14: C:\demo\edtools_demo_data\stagepos_1283\crystal_0001\SMV # Mon Aug 1 21:01:17
 ↪2022

Spgr 1 - Cell 13.64 15.02 25.09 93.07 91.13 114.33 - Vol 1
 ↪4672.25

#	dmax	dmin	ntot	nuniq	compl	i/sig	rmeas	CC(1/2)	ISa	B(ov)	
1	4.35	0.80	583	324	15.0	4.59	13.7	99.0	7.47	6.72	# C:\
↪	0.85	0.80	54	42	12.5	1.96	26.8	91.9			demo\edtools_demo_data\stagepos_0067\crystal_0001\SMV\CORRECT.LP
2	6.39	0.80	229	209	4.8	24.01	22.3	96.3	50.00	4.74	# C:\
↪	0.91	0.85	31	29	4.5	12.16	21.3	0.0			demo\edtools_demo_data\stagepos_0164\crystal_0000\SMV\CORRECT.LP
3	2.05	0.80	400	312	7.5	2.44	20.6	95.2	4.24	6.11	# C:\
↪	0.84	0.80	27	26	4.7	1.79	11.2	0.0			demo\edtools_demo_data\stagepos_0299\crystal_0001\SMV\CORRECT.LP
4	11.09	0.79	3744	2147	42.8	3.44	13.1	99.6	13.90	8.10	# C:\
↪	0.97	0.90	623	336	47.9	1.32	68.9	84.7			demo\edtools_demo_data\stagepos_0325\crystal_0000\SMV\CORRECT.LP
5	6.86	0.80	2161	1081	21.8	4.23	10.5	99.9	33.30	8.74	# C:\
↪	0.97	0.90	342	159	23.3	0.83	130.3	69.6			demo\edtools_demo_data\stagepos_0341\crystal_0000\SMV\CORRECT.LP
6	10.17	0.80	611	443	9.4	3.17	14.9	97.5	5.07	4.64	# C:\
↪	0.85	0.80	56	53	7.0	1.96	73.1	6.4			demo\edtools_demo_data\stagepos_0368\crystal_0001\SMV\CORRECT.LP
7	5.10	0.80	443	330	7.0	3.80	10.7	99.4	8.61	5.62	# C:\
↪											demo\edtools_demo_data\stagepos_0538\crystal_0000\SMV\CORRECT.LP

(continues on next page)

(continued from previous page)

```

- 0.85 0.80 38 36 4.8 1.80 76.5 0.0
8 6.37 0.80 1460 989 9.1 2.88 16.3 97.8 5.24 7.62 # C:\
↪demo\edtools_demo_data\stagepos_0648\crystal_0001\SMV\CORRECT.LP
- 0.85 0.80 166 125 7.3 1.36 62.8 56.0
9 13.11 0.79 2063 1319 22.1 3.46 10.5 99.6 12.09 7.58 # C:\
↪demo\edtools_demo_data\stagepos_0849\crystal_0000\SMV\CORRECT.LP
- 0.89 0.83 326 223 24.5 1.01 53.6 83.4
10 12.33 0.80 479 300 13.8 3.68 13.4 99.6 16.07 9.46 # C:\
↪demo\edtools_demo_data\stagepos_0905\crystal_0000\SMV\CORRECT.LP
- 1.20 1.07 58 35 14.6 4.72 24.9 88.5
11 11.49 0.80 1596 1144 10.7 3.30 12.4 98.5 7.24 7.18 # C:\
↪demo\edtools_demo_data\stagepos_0905\crystal_0001\SMV\CORRECT.LP
- 0.85 0.80 124 121 7.0 0.94 22.6 83.4
12 7.54 0.80 1746 1222 11.3 4.00 13.3 98.6 8.77 5.85 # C:\
↪demo\edtools_demo_data\stagepos_0980\crystal_0000\SMV\CORRECT.LP
- 0.85 0.80 164 146 8.4 1.48 36.5 86.5
13 5.01 0.81 447 328 7.5 4.11 10.9 98.5 6.65 6.67 # C:\
↪demo\edtools_demo_data\stagepos_1014\crystal_0000\SMV\CORRECT.LP
- 0.85 0.80 51 44 6.3 2.11 18.7 92.7
14 6.60 0.80 3124 2149 11.3 3.54 8.4 99.5 12.64 6.94 # C:\
↪demo\edtools_demo_data\stagepos_1283\crystal_0001\SMV\CORRECT.LP
- 0.85 0.80 346 280 9.2 1.24 56.2 84.9

```

Wrote 14 cells to file cells.xlsx

Wrote 14 cells to file cells.yaml

Wrote 8 entries to file filelist.txt (completeness > 10.0%, CC(1/2) > 90.0%)

Most likely lattice types:

```

1 Lattice type `aP` (spgr: 1) was found 9 times (score: 10056)
2 Lattice type `mC` (spgr: 5) was found 4 times (score: 8551)
3 Lattice type `mP` (spgr: 3) was found 1 times (score: 479)


```

** the score corresponds to the total number of indexed reflections.

10.1.3 Unit-cell-based clustering for phase analysis

```
[ ]: !edtools.find_cell cells.yaml -s --cluster --metric lcv
```

```
[3]: from IPython.display import Image
Image('find_cell_step3.png', embed=True)
```



nbsphinx-code-borderwidth=1

Console Output

```

Linkage method = average
Cutoff distance = 0.078
Distance metric = lcv

-----

Cluster #1 (4 items)
  1 [  5.47  14.07  15.30  63.22  87.59  88.58] Vol.: 1050.9
  3 [  5.33  14.99  16.06  64.44  89.16  82.51] Vol.: 1144.9
 10 [  5.05  14.37  14.53  62.13  88.52  89.11] Vol.:  932.0
 13 [  5.30  14.89  15.18  66.79  86.51  86.59] Vol.: 1098.1
---
Mean: [  5.29  14.58  15.27  64.15  87.95  86.70] Vol.: 1056.5
Min:  [  5.05  14.07  14.53  62.13  86.51  82.51] Vol.:  932.0
Max:  [  5.47  14.99  16.06  66.79  89.16  89.11] Vol.: 1144.9

Cluster #2 (3 items)
  2 [  9.52   9.98  12.85  65.60  87.80  85.43] Vol.: 1107.8
  6 [ 10.21  10.36  12.08  85.86  67.02  81.83] Vol.: 1165.3
  7 [ 10.55  10.75  11.75  80.34  66.42  75.73] Vol.: 1179.4
---
Mean: [ 10.09  10.36  12.23  77.27  73.75  81.00] Vol.: 1150.9
Min:  [  9.52   9.98  11.75  65.60  66.42  75.73] Vol.: 1107.8
Max:  [ 10.55  10.75  12.85  85.86  87.80  85.43] Vol.: 1179.4

Cluster #3 (6 items)
  4 [ 14.04  14.39  14.72  76.68  62.79  61.86] Vol.: 2331.3

```

(continues on next page)

(continued from previous page)

```

  5 [ 13.50  14.38  14.63  75.73  64.60  63.07] Vol.: 2283.0
  8 [ 13.89  14.29  17.00  72.43  63.61  63.57] Vol.: 2684.8
  9 [ 14.81  15.07  15.52  62.45  74.78  62.16] Vol.: 2711.1
 11 [ 13.73  14.56  16.03  84.26  68.05  62.57] Vol.: 2629.5
 12 [ 14.43  14.90  15.40  81.24  74.01  61.15] Vol.: 2787.8
---
Mean: [ 14.07  14.60  15.55  75.46  67.97  62.40] Vol.: 2571.3
Min:  [ 13.50  14.29  14.63  62.45  62.79  61.15] Vol.: 2283.0
Max:  [ 14.81  15.07  17.00  84.26  74.78  63.57] Vol.: 2787.8

Wrote cluster 1 to file `cells_cluster_1_4-items.yaml`
Wrote cluster 2 to file `cells_cluster_2_3-items.yaml`
Wrote cluster 3 to file `cells_cluster_3_6-items.yaml`

```

The three resulted clusters 1, 2, 3 correspond to phases ***CTH**, **RTH**, and **IWV**, respectively.

With the averaged primitive unit cell parameters of each cluster, one can use the online tool http://cci.lbl.gov/cctbx/lattice_symmetry.html to find unit cell with higher symmetry with a pre-set tolerance.

We take cluster 3 (phase **IWV**) as an example. The averaged unit cell parameters are: 14.07, 14.6, 15.55, 75.46, 67.97, 62.40

The unit cell parameters with a higher symmetry (space group: *Fmmm* (69)) are: 14.07, 25.8828, 28.9294, 90, 90, 90

The same operation can be done for all the other clusters.

10.1.4 Update the XDS.INP files

This step used `edtools.update_xds` to update the XDS input files with the determined unit cell parameters and space group.

```
[4]: !edtools.update_xds -c 14.07 25.8828 28.9294 90 90 90 -s 69
```

```

16 files named XDS.INP (subdir: None) found.
C:\demo\edtools_demo_data\stagepos_0067\crystal_0001\SMV\XDS.INP
C:\demo\edtools_demo_data\stagepos_0164\crystal_0000\SMV\XDS.INP
C:\demo\edtools_demo_data\stagepos_0290\crystal_0002\SMV\XDS.INP
C:\demo\edtools_demo_data\stagepos_0299\crystal_0001\SMV\XDS.INP
C:\demo\edtools_demo_data\stagepos_0325\crystal_0000\SMV\XDS.INP
C:\demo\edtools_demo_data\stagepos_0341\crystal_0000\SMV\XDS.INP
C:\demo\edtools_demo_data\stagepos_0368\crystal_0001\SMV\XDS.INP
C:\demo\edtools_demo_data\stagepos_0538\crystal_0000\SMV\XDS.INP
C:\demo\edtools_demo_data\stagepos_0648\crystal_0001\SMV\XDS.INP
C:\demo\edtools_demo_data\stagepos_0849\crystal_0000\SMV\XDS.INP
C:\demo\edtools_demo_data\stagepos_0905\crystal_0000\SMV\XDS.INP
C:\demo\edtools_demo_data\stagepos_0905\crystal_0001\SMV\XDS.INP
C:\demo\edtools_demo_data\stagepos_0980\crystal_0000\SMV\XDS.INP
C:\demo\edtools_demo_data\stagepos_1014\crystal_0000\SMV\XDS.INP
C:\demo\edtools_demo_data\stagepos_1261\crystal_0001\SMV\XDS.INP
C:\demo\edtools_demo_data\stagepos_1283\crystal_0001\SMV\XDS.INP
Updated 16 files

```

10.1.5 Refine phases

Rerun **autoindex**, **extract_xds_info** and **find_cell** for the desired phases to be successfully indexed by *XDS*. All the other phases are hopefully excluded in that a phase with different enough unit cell will not be indexed successfully. There are however cases when different phases have similar unit cells, which cannot be told apart during this step.

[5]: !edtools.autoindex

```

!!! ERROR !!! WRONG TYPE OF INPUT FILE SPECIFIED
!!! ERROR !!! WRONG TYPE OF INPUT FILE SPECIFIED
16 files named XDS.INP (subdir: None) found.

  4: C:\demo\edtools_demo_data\stagepos_0325\crystal_0000\SMV # Mon Aug  1 21:03:22
↳2022
Spgr  69 - Cell      13.88   25.44   27.26   90.00   90.00   90.00 - Vol  ↳
↳9625.70

#   dmax  dmin   ntot   nuniq   compl   i/sig   rmeas  CC(1/2)   ISa   B(ov)
4   9.30  0.80   3938   1852    69.0    3.24    20.2   99.4   11.45   8.21
-   0.91  0.85    614    290    74.4    0.86   109.5   81.7

  5: C:\demo\edtools_demo_data\stagepos_0341\crystal_0000\SMV # Mon Aug  1 21:03:24
↳2022
Spgr  69 - Cell      13.52   24.94   27.07   90.00   90.00   90.00 - Vol  ↳
↳9127.70

#   dmax  dmin   ntot   nuniq   compl   i/sig   rmeas  CC(1/2)   ISa   B(ov)
5  10.88  0.80   2203   1029    40.4    3.84    11.7   99.9   27.38   9.78
-   1.07  0.98    299    135    41.8    1.04   107.2   76.5

  8: C:\demo\edtools_demo_data\stagepos_0648\crystal_0001\SMV # Mon Aug  1 21:03:28
↳2022
Spgr  69 - Cell      14.01   25.97   29.04   90.00   90.00   90.00 - Vol  ↳
↳10565.90

#   dmax  dmin   ntot   nuniq   compl   i/sig   rmeas  CC(1/2)   ISa   B(ov)
8   7.14  0.80   1466    781    26.2    2.61    18.2   97.2    4.73    7.15
-   0.84  0.80    142     92    19.9    0.98    62.3   52.7

  9: C:\demo\edtools_demo_data\stagepos_0849\crystal_0000\SMV # Mon Aug  1 21:03:30
↳2022
Spgr  69 - Cell      15.10   26.02   26.72   90.00   90.00   90.00 - Vol  ↳
↳10498.34

#   dmax  dmin   ntot   nuniq   compl   i/sig   rmeas  CC(1/2)   ISa   B(ov)
9   7.24  0.80   1994   1126    38.5    3.27    11.9   99.5   12.91    8.08
-   0.98  0.90    322    166    41.2    1.27    70.2   89.0

 10: C:\demo\edtools_demo_data\stagepos_0905\crystal_0000\SMV -> Error in IDXREF:
↳RETURN CODE IS IER=          0

```

(continues on next page)

(continued from previous page)

```

11: C:\demo\edtools_demo_data\stagepos_0905\crystal_0001\SMV # Mon Aug 1 21:03:32
↪2022
Spgr 69 - Cell      13.83   25.80   28.73   90.00   90.00   90.00 - Vol  ↪
↪10251.27

#   dmax  dmin   ntot   nuniq   compl   i/sig   rmeas  CC(1/2)   ISa   B(ov)
11   7.08  0.80  1591    808    28.2    2.88    17.1   98.1     6.24  7.63
-    0.90  0.85   254    128    30.4    1.17    42.6   95.4

12: C:\demo\edtools_demo_data\stagepos_0980\crystal_0000\SMV # Mon Aug 1 21:03:34
↪2022
Spgr 69 - Cell      14.39   25.16   28.10   90.00   90.00   90.00 - Vol  ↪
↪10173.67

#   dmax  dmin   ntot   nuniq   compl   i/sig   rmeas  CC(1/2)   ISa   B(ov)
12   5.12  0.80  1669    851    30.2    3.75    16.8   98.0     6.26  5.76
-    0.85  0.80   153    109    25.3    1.34    46.1   68.5

15: C:\demo\edtools_demo_data\stagepos_1283\crystal_0001\SMV # Mon Aug 1 21:03:39
↪2022
Spgr 69 - Cell      13.54   25.23   27.30   90.00   90.00   90.00 - Vol  ↪
↪9326.07

#   dmax  dmin   ntot   nuniq   compl   i/sig   rmeas  CC(1/2)   ISa   B(ov)
15   5.97  0.80  1620    563    21.7    6.15     8.4   99.8    11.79  7.17
-    0.85  0.80   187     78    19.1    2.24    45.1   97.9

```

[6]: !edtools.extract_xds_info

```

7 files named CORRECT.LP (subdir: None) found.
1: C:\demo\edtools_demo_data\stagepos_0325\crystal_0000\SMV # Mon Aug 1 21:03:22
↪2022
Spgr 69 - Cell      13.88   25.44   27.26   90.00   90.00   90.00 - Vol  ↪
↪9625.70

2: C:\demo\edtools_demo_data\stagepos_0341\crystal_0000\SMV # Mon Aug 1 21:03:24
↪2022
Spgr 69 - Cell      13.52   24.94   27.07   90.00   90.00   90.00 - Vol  ↪
↪9127.70

3: C:\demo\edtools_demo_data\stagepos_0648\crystal_0001\SMV # Mon Aug 1 21:03:28
↪2022
Spgr 69 - Cell      14.01   25.97   29.04   90.00   90.00   90.00 - Vol  ↪
↪10565.90

4: C:\demo\edtools_demo_data\stagepos_0849\crystal_0000\SMV # Mon Aug 1 21:03:30
↪2022
Spgr 69 - Cell      15.10   26.02   26.72   90.00   90.00   90.00 - Vol  ↪
↪10498.34

```

(continues on next page)

(continued from previous page)

```

5: C:\demo\edtools_demo_data\stagepos_0905\crystal_0001\SMV # Mon Aug 1 21:03:32
↪2022
Spgr 69 - Cell      13.83   25.80   28.73   90.00   90.00   90.00 - Vol  ↪
↪10251.27

```

```

6: C:\demo\edtools_demo_data\stagepos_0980\crystal_0000\SMV # Mon Aug 1 21:03:34
↪2022
Spgr 69 - Cell      14.39   25.16   28.10   90.00   90.00   90.00 - Vol  ↪
↪10173.67

```

```

7: C:\demo\edtools_demo_data\stagepos_1283\crystal_0001\SMV # Mon Aug 1 21:03:39
↪2022
Spgr 69 - Cell      13.54   25.23   27.30   90.00   90.00   90.00 - Vol  ↪
↪9326.07

```

#	dmax	dmin	ntot	nuniq	compl	i/sig	rmeas	CC(1/2)	ISa	B(ov)	
1	9.30	0.80	3938	1852	69.0	3.24	20.2	99.4	11.45	8.21	# C:\
↪	0.91	0.85	614	290	74.4	0.86	109.5	81.7			demo\edtools_demo_data\stagepos_0325\crystal_0000\SMV\CORRECT.LP
2	10.88	0.80	2203	1029	40.4	3.84	11.7	99.9	27.38	9.78	# C:\
↪	1.07	0.98	299	135	41.8	1.04	107.2	76.5			demo\edtools_demo_data\stagepos_0341\crystal_0000\SMV\CORRECT.LP
3	7.14	0.80	1466	781	26.2	2.61	18.2	97.2	4.73	7.15	# C:\
↪	0.84	0.80	142	92	19.9	0.98	62.3	52.7			demo\edtools_demo_data\stagepos_0648\crystal_0001\SMV\CORRECT.LP
4	7.24	0.80	1994	1126	38.5	3.27	11.9	99.5	12.91	8.08	# C:\
↪	0.98	0.90	322	166	41.2	1.27	70.2	89.0			demo\edtools_demo_data\stagepos_0849\crystal_0000\SMV\CORRECT.LP
5	7.08	0.80	1591	808	28.2	2.88	17.1	98.1	6.24	7.63	# C:\
↪	0.90	0.85	254	128	30.4	1.17	42.6	95.4			demo\edtools_demo_data\stagepos_0905\crystal_0001\SMV\CORRECT.LP
6	5.12	0.80	1669	851	30.2	3.75	16.8	98.0	6.26	5.76	# C:\
↪	0.85	0.80	153	109	25.3	1.34	46.1	68.5			demo\edtools_demo_data\stagepos_0980\crystal_0000\SMV\CORRECT.LP
7	5.97	0.80	1620	563	21.7	6.15	8.4	99.8	11.79	7.17	# C:\
↪	0.85	0.80	187	78	19.1	2.24	45.1	97.9			demo\edtools_demo_data\stagepos_1283\crystal_0001\SMV\CORRECT.LP

Wrote 7 cells to file cells.xlsx

Wrote 7 cells to file cells.yaml

Wrote 7 entries to file filelist.txt (completeness > 10.0%, CC(1/2) > 90.0%)

Most likely lattice types:


(continues on next page)

(continued from previous page)

```
1 Lattice type `oF` (spgr: 22) was found 7 times (score: 14481)
** the score corresponds to the total number of indexed reflections.
```

```
[ ]: !edtools.find_cell cells.yaml --cluster --metric lcv
```

```
[7]: Image('find_cell_step5.png', embed=True)
```



nbsphinx-code-borderw

Console Output

```
Linkage method = average
Cutoff distance = 0.0551
Distance metric = lcv
```

Cluster #1 (7 items)

1	[13.97	25.49	27.12	90.00	90.00	90.00]	Vol.: 9657.9
2	[13.53	25.01	27.18	90.00	90.00	90.00]	Vol.: 9195.6
3	[14.03	26.02	29.55	90.00	90.00	90.00]	Vol.: 10790.3
4	[14.94	26.14	26.94	90.00	90.00	90.00]	Vol.: 10522.3
5	[13.85	25.79	29.03	90.00	90.00	90.00]	Vol.: 10364.0
6	[14.52	24.95	28.11	90.00	90.00	90.00]	Vol.: 10184.6
7	[13.53	25.13	27.15	90.00	90.00	90.00]	Vol.: 9233.7

Mean:	[14.05	25.50	27.87	90.00	90.00	90.00]	Vol.: 9992.6
Min:	[13.53	24.95	26.94	90.00	90.00	90.00]	Vol.: 9195.6
Max:	[14.94	26.14	29.55	90.00	90.00	90.00]	Vol.: 10790.3

(continues on next page)

(continued from previous page)

```
Wrote cluster 1 to file `cells_cluster_1_7-items.yaml`
```

10.1.6 Generate the input file for *XSCALE*

This command generates the desired unit cell cluster for *XSCALE*.

```
[8]: !edtools.make_xscale cells_cluster_1_7-items.yaml -c 14.05 25.50 27.87 90.00 90.00 90.00
↵-s 69

Loaded 7 cells
Lowest possible symmetry for 69 (oF): 22

Using:
  SPACE_GROUP_NUMBER= 69
  UNIT_CELL_CONSTANTS= 14.050 25.500 27.870 90.000 90.000 90.000

Wrote file XSCALE.INP
Wrote file XDSCONV.INP
```

10.1.7 Run *XSCALE*

XSCALE calculates the correlation coefficients between different datasets.

```
[9]: !wsl xscale

***** XSCALE ***** (VERSION Jan 10, 2022 BUILT=20220220) 1-Aug-2022
Author: Wolfgang Kabsch
Copy licensed until 31-Mar-2023 to
  academic users for non-commercial applications
No redistribution.

*****
CONTROL CARDS
*****

SNRC= 2
SAVE_CORRECTION_IMAGES= FALSE
SPACE_GROUP_NUMBER= 69
UNIT_CELL_CONSTANTS= 14.050 25.500 27.870 90.000 90.000 90.000

OUTPUT_FILE= MERGED.HKL

INPUT_FILE= edtools_demo_data/stagepos_0325/crystal_0000/SMV/XDS_ASCII.HKL
INCLUDE_RESOLUTION_RANGE= 20 0.8

INPUT_FILE= edtools_demo_data/stagepos_0341/crystal_0000/SMV/XDS_ASCII.HKL
INCLUDE_RESOLUTION_RANGE= 20 0.8

INPUT_FILE= edtools_demo_data/stagepos_0648/crystal_0001/SMV/XDS_ASCII.HKL
```

(continues on next page)

(continued from previous page)

```
INCLUDE_RESOLUTION_RANGE= 20 0.8
```

```
INPUT_FILE= edtools_demo_data/stagepos_0849/crystal_0000/SMV/XDS_ASCII.HKL
INCLUDE_RESOLUTION_RANGE= 20 0.8
```

```
INPUT_FILE= edtools_demo_data/stagepos_0905/crystal_0001/SMV/XDS_ASCII.HKL
INCLUDE_RESOLUTION_RANGE= 20 0.8
```

```
INPUT_FILE= edtools_demo_data/stagepos_0980/crystal_0000/SMV/XDS_ASCII.HKL
INCLUDE_RESOLUTION_RANGE= 20 0.8
```

```
INPUT_FILE= edtools_demo_data/stagepos_1283/crystal_0001/SMV/XDS_ASCII.HKL
INCLUDE_RESOLUTION_RANGE= 20 0.8
```

THE DATA COLLECTION STATISTICS REPORTED BELOW ASSUMES:

SPACE_GROUP_NUMBER= 69

UNIT_CELL_CONSTANTS= 14.05 25.50 27.87 90.000 90.000 90.000

ALL DATA SETS WILL BE SCALED TO edtools_demo_data/stagepos_0325/crystal_0000/SMV/XDS_
↪ASCII.HKL

```
*****
                        READING INPUT REFLECTION DATA FILES
*****
```

DATA SET#	MEAN INTENSITY	REFLECTIONS ACCEPTED	REFLECTIONS REJECTED	INPUT FILE NAME
1	0.3010E+02	3938	0	edtools_demo_data/stagepos_0325/crystal_0000/SMV/XDS_↪ASCII.HKL
2	0.1368E+02	2205	0	edtools_demo_data/stagepos_0341/crystal_0000/SMV/XDS_↪ASCII.HKL
3	0.9168E+02	1453	0	edtools_demo_data/stagepos_0648/crystal_0001/SMV/XDS_↪ASCII.HKL
4	0.4279E+02	1931	0	edtools_demo_data/stagepos_0849/crystal_0000/SMV/XDS_↪ASCII.HKL
5	0.8542E+02	1590	0	edtools_demo_data/stagepos_0905/crystal_0001/SMV/XDS_↪ASCII.HKL
6	0.1676E+03	1662	0	edtools_demo_data/stagepos_0980/crystal_0000/SMV/XDS_↪ASCII.HKL
7	0.1915E+03	1620	0	edtools_demo_data/stagepos_1283/crystal_0001/SMV/XDS_↪ASCII.HKL

```
*****
                        OVERALL SCALING AND CRYSTAL DISORDER CORRECTION
*****
```

CORRELATIONS BETWEEN INPUT DATA SETS AFTER CORRECTIONS

(continues on next page)

(continued from previous page)

DATA SETS #i	NUMBER OF COMMON #j	REFLECTIONS	CORRELATION BETWEEN i,j	RATIO OF COMMON INTENSITIES (i/j)	B-FACTOR BETWEEN i,j
1	2	119	0.936	3.1896	0.1394
1	3	87	0.980	0.5337	-0.4854
2	3	164	0.945	0.1982	-0.5975
1	4	216	0.925	1.0428	-0.6421
2	4	116	0.972	0.3076	-0.4582
3	4	131	0.894	1.6716	0.0043
1	5	80	0.959	0.3779	0.1598
2	5	147	0.970	0.1872	-0.3822
3	5	218	0.988	0.9850	0.0336
4	5	96	0.928	0.6114	0.1442
1	6	206	0.955	0.1970	-0.4620
2	6	81	0.949	0.0917	-0.7741
3	6	91	0.934	0.6363	-0.3722
4	6	106	0.927	0.2401	-0.0297
5	6	81	0.866	0.5334	-0.1957
1	7	35	0.965	0.4893	-1.1586
2	7	131	0.981	0.1113	-1.0633
3	7	158	0.984	0.5449	-0.3181
4	7	45	0.833	0.4069	-0.1510
5	7	196	0.987	0.5649	-0.4493
6	7	67	0.846	1.6099	-0.5928

K*EXP(B*SS) = Factor applied to intensities
 SS = (2sin(theta)/lambda)^2

K	B	DATA SET NAME
1.000E+00	0.000	edtools_demo_data/stagepos_0325/crystal_0000/SMV/XDS_ASCII.HKL
2.961E+00	0.170	edtools_demo_data/stagepos_0341/crystal_0000/SMV/XDS_ASCII.HKL
5.374E-01	-0.365	edtools_demo_data/stagepos_0648/crystal_0001/SMV/XDS_ASCII.HKL
9.426E-01	-0.465	edtools_demo_data/stagepos_0849/crystal_0000/SMV/XDS_ASCII.HKL
5.063E-01	-0.243	edtools_demo_data/stagepos_0905/crystal_0001/SMV/XDS_ASCII.HKL
2.304E-01	-0.530	edtools_demo_data/stagepos_0980/crystal_0000/SMV/XDS_ASCII.HKL
3.491E-01	-0.812	edtools_demo_data/stagepos_1283/crystal_0001/SMV/XDS_ASCII.HKL

 CORRECTION PARAMETERS FOR THE STANDARD ERROR OF REFLECTION INTENSITIES

The variance $v_0(I)$ of the intensity I obtained from counting statistics is replaced by $v(I)=a*(v_0(I)+b*I^2)$. The model parameters a , b are chosen to minimize the discrepancies between $v(I)$ and the variance estimated from sample statistics of symmetry related reflections. This model implicates an asymptotic limit $ISa=1/\text{SQRT}(a*b)$ for the highest $I/\text{Sigma}(I)$ that the experimental setup can produce (Diederichs (2010) Acta Cryst D66, 733-740). Often the value of ISa is reduced from the initial value ISa_0 due to systematic errors showing up by comparison with other data sets in the scaling procedure. ($ISa=ISa_0=-1$ if v_0 is unknown for a data set.)

(continues on next page)

(continued from previous page)

a	b	ISa	ISa0	INPUT DATA SET
3.014E+00	1.258E-02	5.14	11.45	edtools_demo_data/stagepos_0325/crystal_0000/SMV/ ↪XDS_ASCII.HKL
2.201E+00	3.743E-03	11.02	27.38	edtools_demo_data/stagepos_0341/crystal_0000/SMV/ ↪XDS_ASCII.HKL
8.809E+00	2.191E-02	2.28	4.73	edtools_demo_data/stagepos_0648/crystal_0001/SMV/ ↪XDS_ASCII.HKL
6.242E+00	1.032E-02	3.94	12.91	edtools_demo_data/stagepos_0849/crystal_0000/SMV/ ↪XDS_ASCII.HKL
7.668E+00	1.817E-02	2.68	6.25	edtools_demo_data/stagepos_0905/crystal_0001/SMV/ ↪XDS_ASCII.HKL
1.379E+01	1.128E-02	2.53	6.26	edtools_demo_data/stagepos_0980/crystal_0000/SMV/ ↪XDS_ASCII.HKL
7.838E-01	1.921E-01	2.58	11.79	edtools_demo_data/stagepos_1283/crystal_0001/SMV/ ↪XDS_ASCII.HKL

FACTOR TO PLACE ALL DATA SETS TO AN APPROXIMATE ABSOLUTE SCALE 0.143057E+03
(ASSUMING A PROTEIN WITH 50% SOLVENT)

STATISTICS OF SCALED OUTPUT DATA SET : MERGED.HKL

FILE TYPE: XDS_ASCII MERGE=FALSE FRIEDEL'S_LAW=TRUE

13 OUT OF 14399 REFLECTIONS REJECTED
14386 REFLECTIONS ON OUTPUT FILE

DEFINITIONS:

R-FACTOR

observed = $(\text{SUM}(\text{ABS}(I(h,i) - I(h)))) / (\text{SUM}(I(h,i)))$

expected = expected R-FACTOR derived from $\text{Sigma}(I)$

COMPARED = number of reflections used for calculating R-FACTOR

I/SIGMA = mean of intensity/ $\text{Sigma}(I)$ of unique reflections
(after merging symmetry-related observations)

$\text{Sigma}(I)$ = standard deviation of reflection intensity I
estimated from sample statistics

R-meas = redundancy independent R-factor (intensities)
Diederichs & Karplus (1997), Nature Struct. Biol. 4, 269-275.

CC(1/2) = percentage of correlation between intensities from
random half-datasets. Correlation significant at
the 0.1% level is marked by an asterisk.
Karplus & Diederichs (2012), Science 336, 1030-33

Anomal
Corr = percentage of correlation between random half-sets
of anomalous intensity differences. Correlation
significant at the 0.1% level is marked.

(continues on next page)

(continued from previous page)

SigAno = mean anomalous difference in units of its estimated standard deviation ($|F(+)-F(-)|/\text{Sigma}$). $F(+)$, $F(-)$ are structure factor estimates obtained from the merged intensity observations in each parity class.

Nano = Number of unique reflections used to calculate Anomal_Corr & SigAno. At least two observations for each (+ and -) parity are required.

cpu time used by XSCALE 0.2 sec
 elapsed wall-clock time 0.2 sec


10.1.8 Intensity-based clustering

Run intensity-based clustering to further filter out datasets with low correlation (to remove poor quality datasets), or from a different phase that with similar enough unit cell. Cut-off on the dendrogram is selected manually. A number below 0.4 can be a good starting choice.

In the end, integration results from datasets corresponding to different clusters are automatically copied to different folders after running clustering. The merged intensities in file `shelx.hkl` can be used for structure determination.

```
[ ]: !edtools.cluster
```

```
[10]: Image('intensity_cluster.png', embed=True)
```



nbsphinx-code-border [\[10\]](#)

Console Output

```
Running XSCALE on cluster 1
Running XSCALE on cluster 2
```

```
Clustering results
```

```
Cutoff distance: 0.252
Equivalent CC(I): 0.968
Method: average
```

#	N_clust	CC(1/2)	N_obs	N_uniq	N_poss	Compl.	N_comp	R_meas	d_min
↔ i/sigma		Lauegr.	prob.	conf.	idx				
2**	2	99.8*	4111	1546	2789	55.4	3723	0.143*	0.80
↔	3.27								
1***	4	97.3*	8599	2496	2782	89.7*	8220	0.270*	0.80
↔	2.85								

(Sorted by 'Completeness')

```
Cluster 1: [1, 3, 5, 7]
Cluster 2: [2, 4]
```

10.2 Instruction for using on your own data

- Install **edtools** and all software dependencies on your system
- Put all your 3D ED datasets in one folder. All the 3D ED datasets are expected to be in some *XDS* readable image format, e.g. *SMV*. A correctly configured *XDS.INP* file is also expected for each dataset.
- Open Windows command prompt from the root directory which contains all the datasets
- Follow the demo

CHAPTER
ELEVEN

LINKS

INDICES AND TABLES

- `genindex`
- `modindex`
- `search`

PYTHON MODULE INDEX

e

- edtools, 39
- edtools.autoindex, 19
- edtools.cif_tools, 20
- edtools.cluster, 20
- edtools.extract_xds_info, 21
- edtools.find_beam_center, 23
- edtools.find_cell, 24
- edtools.find_rotation_axis, 26
- edtools.make_xscale, 27
- edtools.period_table, 28
- edtools.reflection_tool, 29
- edtools.run_pointless, 32
- edtools.scattering_factor, 32
- edtools.update_xds, 36
- edtools.utils, 36
- edtools.widgets, 37
- edtools.wsl, 39

Symbols

- `__init__()` (*edtools.extract_xds_info.xds_parser* method), 22
 - `__init__()` (*edtools.period_table.App* method), 28
 - `__init__()` (*edtools.reflection_tool.GroupReflectionsGUI* method), 30
 - `__init__()` (*edtools.scattering_factor.ScatteringFactorGUI* method), 35
 - `__init__()` (*edtools.widgets Hoverbox* method), 38
 - `__init__()` (*edtools.widgets.Spinbox* method), 38
- ## A
- `a0` (*edtools.scattering_factor.FitResult_4p* property), 33
 - `a0` (*edtools.scattering_factor.FitResult_5p* property), 33
 - `a1` (*edtools.scattering_factor.FitResult_4p* property), 33
 - `a1` (*edtools.scattering_factor.FitResult_5p* property), 34
 - `a2` (*edtools.scattering_factor.FitResult_4p* property), 33
 - `a2` (*edtools.scattering_factor.FitResult_5p* property), 34
 - `a3` (*edtools.scattering_factor.FitResult_4p* property), 33
 - `a3` (*edtools.scattering_factor.FitResult_5p* property), 34
 - `a4` (*edtools.scattering_factor.FitResult_5p* property), 34
 - `add_instrument()` (in module *edtools.cif_tools*), 20
 - `App` (class in *edtools.period_table*), 28
- ## B
- `b0` (*edtools.scattering_factor.FitResult_4p* property), 33
 - `b0` (*edtools.scattering_factor.FitResult_5p* property), 34
 - `b1` (*edtools.scattering_factor.FitResult_4p* property), 33
 - `b1` (*edtools.scattering_factor.FitResult_5p* property), 34
 - `b2` (*edtools.scattering_factor.FitResult_4p* property), 33
 - `b2` (*edtools.scattering_factor.FitResult_5p* property), 34
 - `b3` (*edtools.scattering_factor.FitResult_4p* property), 33
 - `b3` (*edtools.scattering_factor.FitResult_5p* property), 34
 - `b4` (*edtools.scattering_factor.FitResult_5p* property), 34
- ## C
- `c` (*edtools.scattering_factor.FitResult_4p* property), 33
 - `c` (*edtools.scattering_factor.FitResult_5p* property), 34
 - `calc_structure_factors()` (*edtools.reflection_tool.GroupReflectionsGUI* method), 31
 - `cell_as_dict()` (*edtools.extract_xds_info.xds_parser* method), 22
 - `cell_info()` (*edtools.extract_xds_info.xds_parser* method), 22
 - `cells_to_cellparm()` (in module *edtools.extract_xds_info*), 21
 - `cells_to_excel()` (in module *edtools.extract_xds_info*), 21
 - `cells_to_yaml()` (in module *edtools.extract_xds_info*), 22
 - `cells_to_yaml_xparm()` (in module *edtools.extract_xds_info*), 22
 - `check_for_pointless()` (in module *edtools.cluster*), 20
 - `check_I_frame_seq()` (*edtools.reflection_tool.GroupReflectionsGUI* method), 31
 - `clean_params()` (in module *edtools.cluster*), 21
 - `clear_files()` (in module *edtools.autoindex*), 19
 - `clear_plot()` (*edtools.scattering_factor.ScatteringFactorGUI* method), 35
 - `cluster_cell()` (in module *edtools.find_cell*), 24
 - `connect()` (in module *edtools.autoindex*), 19
 - `corr_prec()` (*edtools.reflection_tool.GroupReflectionsGUI* method), 31
 - `cylinder_histo()` (in module *edtools.find_rotation_axis*), 26
- ## D
- `d_calculator()` (in module *edtools.find_cell*), 24
 - `delete_instrument()` (in module *edtools.cif_tools*), 20
 - `distance_from_dendrogram()` (in module *edtools.cluster*), 21
 - `distance_from_dendrogram()` (in module *edtools.find_cell*), 24
 - `draw()` (*edtools.scattering_factor.ScatteringFactorGUI* method), 35
- ## E
- `edtools` module, 39
 - `edtools.autoindex`

- module, 19
 - edtools.cif_tools
 - module, 20
 - edtools.cluster
 - module, 20
 - edtools.extract_xds_info
 - module, 21
 - edtools.find_beam_center
 - module, 23
 - edtools.find_cell
 - module, 24
 - edtools.find_rotation_axis
 - module, 26
 - edtools.make_xscale
 - module, 27
 - edtools.period_table
 - module, 28
 - edtools.reflection_tool
 - module, 29
 - edtools.run_pointless
 - module, 32
 - edtools.scattering_factor
 - module, 32
 - edtools.update_xds
 - module, 36
 - edtools.utils
 - module, 36
 - edtools.widgets
 - module, 37
 - edtools.wsl
 - module, 39
 - enter() (*edtools.widgets Hoverbox* method), 38
 - evaluate_symmetry() (in module *edtools.extract_xds_info*), 22
 - exti_corr() (*edtools.reflection_tool.GroupReflectionsGUI* method), 31
 - exti_corr_der() (*edtools.reflection_tool.GroupReflectionsGUI* method), 31
- ## F
- f_calc_structure_factors() (*edtools.reflection_tool.GroupReflectionsGUI* method), 31
 - find_beam_center() (in module *edtools.find_beam_center*), 23
 - find_cell() (in module *edtools.find_cell*), 25
 - find_peak_max() (in module *edtools.find_beam_center*), 23
 - fit_4_param() (*edtools.scattering_factor.ScatteringFactorGUI* method), 35
 - fit_4_param() (in module *edtools.scattering_factor*), 36
 - fit_5_param() (*edtools.scattering_factor.ScatteringFactorGUI* method), 35
 - fit_5_param() (in module *edtools.scattering_factor*), 36
 - FitResult_4p (class in *edtools.scattering_factor*), 32
 - FitResult_5p (class in *edtools.scattering_factor*), 33
 - func_4p() (*edtools.scattering_factor.ScatteringFactorGUI* method), 35
 - func_4p_charge() (*edtools.scattering_factor.ScatteringFactorGUI* method), 35
 - func_5p() (*edtools.scattering_factor.ScatteringFactorGUI* method), 35
 - func_5p_charge() (*edtools.scattering_factor.ScatteringFactorGUI* method), 36
- ## G
- gather_xds_ascii() (in module *edtools.extract_xds_info*), 22
 - gen_fcalc() (*edtools.reflection_tool.GroupReflectionsGUI* method), 31
 - get_clusters() (in module *edtools.cluster*), 21
 - get_clusters() (in module *edtools.find_cell*), 25
 - get_condensed_distance_matrix() (in module *edtools.cluster*), 21
 - get_parameters() (*edtools.scattering_factor.ScatteringFactorGUI* method), 36
 - get_xds_ascii_names() (in module *edtools.make_xscale*), 27
 - group_df() (*edtools.reflection_tool.GroupReflectionsGUI* method), 31
 - GroupReflectionsGUI (class in *edtools.reflection_tool*), 29
- ## H
- hidetip() (*edtools.widgets Hoverbox* method), 38
 - Hoverbox (class in *edtools.widgets*), 37
- ## I
- info() (*edtools.period_table.App* method), 28
 - info_header() (*edtools.extract_xds_info.xds_parser* static method), 22
 - init_vars() (*edtools.reflection_tool.GroupReflectionsGUI* method), 31
 - init_vars() (*edtools.scattering_factor.ScatteringFactorGUI* method), 36
 - integration_info() (*edtools.extract_xds_info.xds_parser* method), 23

L

`lattice_to_space_group()` (in module `edtools.extract_xds_info`), 22

`leave()` (`edtools.widgets.Hoverbox` method), 38

`list_instrument()` (in module `edtools.cif_tools`), 20

`load_lib()` (`edtools.scattering_factor.ScatteringFactorGUI` method), 36

`load_spot_xds()` (in module `edtools.find_rotation_axis`), 26

`lorentz_corr()` (`edtools.reflection_tool.GroupReflectionsGUI` method), 31

M

`main()` (in module `edtools.autoindex`), 19

`main()` (in module `edtools.cif_tools`), 20

`main()` (in module `edtools.cluster`), 21

`main()` (in module `edtools.extract_xds_info`), 22

`main()` (in module `edtools.find_beam_center`), 23

`main()` (in module `edtools.find_cell`), 25

`main()` (in module `edtools.find_rotation_axis`), 26

`main()` (in module `edtools.make_xscale`), 27

`main()` (in module `edtools.period_table`), 28

`main()` (in module `edtools.reflection_tool`), 32

`main()` (in module `edtools.run_pointless`), 32

`main()` (in module `edtools.scattering_factor`), 36

`main()` (in module `edtools.update_xds`), 36

`make()` (in module `edtools.find_rotation_axis`), 26

`make_2d_rotmat()` (in module `edtools.find_rotation_axis`), 27

module

`edtools`, 39

`edtools.autoindex`, 19

`edtools.cif_tools`, 20

`edtools.cluster`, 20

`edtools.extract_xds_info`, 21

`edtools.find_beam_center`, 23

`edtools.find_cell`, 24

`edtools.find_rotation_axis`, 26

`edtools.make_xscale`, 27

`edtools.period_table`, 28

`edtools.reflection_tool`, 29

`edtools.run_pointless`, 32

`edtools.scattering_factor`, 32

`edtools.update_xds`, 36

`edtools.utils`, 36

`edtools.widgets`, 37

`edtools.wsl`, 39

N

`name()` (`edtools.period_table.App` method), 28

O

`open_file()` (`edtools.reflection_tool.GroupReflectionsGUI` method), 31

`optimize()` (in module `edtools.find_rotation_axis`), 27

`output()` (`edtools.scattering_factor.ScatteringFactorGUI` method), 36

P

`parse()` (`edtools.extract_xds_info.xds_parser` method), 23

`parse_args_for_fns()` (in module `edtools.utils`), 37

`parse_cellparm()` (in module `edtools.find_cell`), 25

`parse_xds()` (in module `edtools.autoindex`), 19

`parse_xds_ascii()` (in module `edtools.make_xscale`), 28

`parse_xds_inp()` (in module `edtools.find_rotation_axis`), 27

`parse_xparm_for_uc()` (in module `edtools.extract_xds_info`), 22

`parse_xscale_lp()` (in module `edtools.cluster`), 21

`parse_xscale_lp_initial()` (in module `edtools.cluster`), 21

`plot()` (`edtools.scattering_factor.ScatteringFactorGUI` method), 36

`plot_histo()` (in module `edtools.find_rotation_axis`), 27

`point_to_rotation_axis()` (`edtools.reflection_tool.GroupReflectionsGUI` method), 31

`print_filename()` (`edtools.extract_xds_info.xds_parser` method), 23

`put_in_order()` (in module `edtools.find_cell`), 25

R

`random_sample()` (in module `edtools.find_rotation_axis`), 27

`read()` (`edtools.scattering_factor.ScatteringFactorGUI` method), 36

`read_adsc()` (in module `edtools.find_beam_center`), 23

`read_cif()` (`edtools.reflection_tool.GroupReflectionsGUI` method), 31

`readheader()` (in module `edtools.find_beam_center`), 23

`remove_reflection()` (`edtools.reflection_tool.GroupReflectionsGUI` method), 31

`rotation_axis_to_xyz()` (in module `edtools.find_rotation_axis`), 27

`rotation_matrix()` (in module `edtools.find_rotation_axis`), 27

`run_pointless()` (in module `edtools.cluster`), 21

`run_xscale()` (in module `edtools.cluster`), 21

S

`save_file()` (`edtools.reflection_tool.GroupReflectionsGUI` method), 31

- save_grouped() (*edtools.reflection_tool.GroupReflectionsGUI* method), 31
- scaling_factor() (*edtools.reflection_tool.GroupReflectionsGUI* method), 31
- scaling_func() (*edtools.reflection_tool.GroupReflectionsGUI* method), 32
- ScatteringFactorGUI (class in *edtools.scattering_factor*), 34
- schedule() (*edtools.widgets.Hoverbox* method), 38
- select_element() (*edtools.scattering_factor.ScatteringFactorGUI* method), 36
- set() (*edtools.widgets.Spinbox* method), 38
- showtip() (*edtools.widgets.Hoverbox* method), 38
- space_group (*edtools.extract_xds_info.xds_parser* property), 23
- space_group_lib() (in module *edtools.utils*), 37
- Spinbox (class in *edtools.widgets*), 38
- split_grouped() (*edtools.reflection_tool.GroupReflectionsGUI* method), 32
- swap_needed() (in module *edtools.find_beam_center*), 24
- ## T
- to_radian() (in module *edtools.find_cell*), 25
- to_sin() (in module *edtools.find_cell*), 25
- transform_integrated() (*edtools.reflection_tool.GroupReflectionsGUI* method), 32
- translate_image() (in module *edtools.find_beam_center*), 24
- ## U
- unit_cell (*edtools.extract_xds_info.xds_parser* property), 23
- unit_cell_lcv_distance() (in module *edtools.find_cell*), 25
- unsubscribe() (*edtools.widgets.Hoverbox* method), 38
- update_cif() (in module *edtools.cif_tools*), 20
- update_xds() (in module *edtools.update_xds*), 36
- ## V
- validate() (*edtools.reflection_tool.GroupReflectionsGUI* method), 32
- validate_range() (*edtools.reflection_tool.GroupReflectionsGUI* method), 32
- volume (*edtools.extract_xds_info.xds_parser* property), 23
- volume() (in module *edtools.utils*), 37
- volume_difference() (in module *edtools.find_cell*), 25
- weighted_average() (in module *edtools.find_cell*), 25
- write_adsc() (in module *edtools.find_beam_center*), 24
- write_xdsconv_inp() (in module *edtools.make_xscale*), 28
- write_xscale_inp() (in module *edtools.make_xscale*), 28
- ## X
- xds_index() (in module *edtools.autoindex*), 19
- xds_parser (class in *edtools.extract_xds_info*), 22
- xyz2cyl() (in module *edtools.find_rotation_axis*), 27